



# Product information

## Drill Bath PLUS

2500 ml    SUN410K

10000 ml    SUN410N



## Characteristics

- Ready-for-use universal solution for the easy preparation of drills and rotating instruments.
- Due to the mild pH value and suitable substances, the sharpness of the cutting edges of drills is excellently preserved.
- Good corrosion protection, high material compatibility, and pleasant smell.

## Product type

Highly effective and cleansing, ready-for-use combination cleaning product with a comprehensive range of efficacy against bacteria (including TB), yeast and enveloped viruses. The product provides effective corrosion protection.

## Areas of application

For the cleaning and prophylactic disinfection of drills, cutters, diamond cutters and other rotating instruments. Not suitable for the terminal disinfection of semi-critical and critical medical devices!

## Application

Always use the PLUS drill bath undiluted! After use, put the instruments straight into the ready-for-use solution, all the parts to be disinfected have to be fully covered. After disinfection the instruments should be rinsed, ideally with demineralised water, and dried thoroughly. Do not mix with other disinfectants or cleaners, not suitable for aluminium. Depending on the degree of contamination, the working solution may need to be replaced (RKI recommendation: at least once per working day).

The working solution can also be disposed of in the sewage water system.

It is recommended to wear safety gloves and eye protection.

## Chemical-physical data

Density: 1.01g/cm<sup>3</sup>

pH value: 9.8

## Available containers

- Bottles and canisters

## Application times

Bacteria (incl. MRSA) as per EN 13727, EN 14561 and VAH, dirty conditions:

Conc.: 5 min.

Levurocidal (Candida albicans) as per EN 13624, EN 14562 and VAH, high contamination:

Conc.: 5 min.

Tuberculocidal (*M. terrae*) as per EN 14348, EN 14563, clean conditions:

Conc.: 5 min.

Limited virucidal effect (enveloped viruses e.g. HBV, HCV and HIV) as per RKI/DVV (dirty conditions):

Conc.: 5 min. (dirty conditions)

## Extra application times

Bacteria as per EN 13727 (dirty conditions):

100% 1 min.

Candida albicans as per EN 13624 (dirty conditions):

100% 1 min.

Bacteria as per EN 14561 (dirty conditions):

100% 5 min.

Candida albicans as per EN 14562 (dirty conditions):

100% 1 min.

Bacteria and candida albicans as per VAH (high contamination):

100% 5 min.

TbB (*M. terrae*) as per EN 14348, EN 14563 (clean conditions):

100% 5 min.

"Limited virucidal effect" (enveloped viruses as per RKI/DVV such as HBV, HCV and HIV), dirty conditions:

12% 5 min.

## Safety guidelines

H412 Harmful to aquatic organisms, with long-term effect. Avoid releasing into the environment. Dispose of the contents/container in accordance with local/regional/national/international regulations.

## Contact

PRISMAN Pharma International AG

Am Stalden 16

4622 Egerkingen

Switzerland

Tel.: +41 (0)62 206 20 60

Fax: +41 (0)62 206 20 69

[www.innocid.com](http://www.innocid.com)

[info@prismanpharma.com](mailto:info@prismanpharma.com)

## Manufacturer

PRISMAN GmbH

Otto-Hahn-Ring 6-18

64653 Lorsch

Germany