



# Product information

## Surface Disinfection PLUS Konz.

<b>1000 ml</b>	<b>SUN370I</b>
<b>2500 ml</b>	<b>SUN370K</b>
<b>5000 ml</b>	<b>SUN370L</b>
<b>10000 ml</b>	<b>SUN370N</b>



- Concentrate for the regular cleaning and disinfection of floors, other surfaces and medical products
- Suitable for dispenser equipment
- Extremely wide spectrum of efficacy for daily demands
- Strong grease-removing effect
- Lowest application concentrations (super concentrate)
- Very mild on all surfaces
- Leaves no sticky residue and does not foam
- Licensed as a medical product and biocide

## Product type

Water-based disinfectant concentrate for the cleaning and disinfection of surfaces. In the event of low application concentrations, there is an extensive spectrum of action against bacteria, yeast and viruses.

## Areas of application

### As per the Biocidal Products Directive:

Aldehyde-free surface concentrate for the disinfection and cleaning of surfaces.

### As per the Medical Devices Law (MPG):

Surface disinfectant for the prophylactic surface disinfection of medical products. Not suitable for the terminal disinfection of invasive medical products.

## Application

Before use, the disinfectant is diluted with water. Do not mix with cleaning agents or other disinfectants! Before using the product on large areas of sensitive materials like artificial leather, first test the material tolerance on a small, hidden area. To prevent streaks forming, the treated areas should be wiped dry at regular intervals with a damp (with water) disposable cloth.

## Application times

Bacteria and candida albicans (EN 13727, EN 13697, EN 13624, VAH, and EN 16615, dirty conditions):

## Composition

100g contains:

10.8g didecyldimethylammonium chloride

6.0g N-(3-aminopropyl)-N-dodecylpropane-1.3-diamin.

## Chemical-physical data

Density: 1.01g/m<sup>3</sup>

pH-value (at 20° C): >11

## Supplementary products

- Disinfection wipes
- Alcoholic spray disinfectants
- Accessories (canister pump, canister tap, canister key etc.)

0.5% 15 min.

TbB (*m. terrae*) as per EN 14348:

2.0% 15 min. clean conditions

0.5% 60 min. clean conditions

1.0% 60 min, dirty conditions

"Limited virucidal effect" (enveloped viruses as per RKI/DVV such as HBV, HCV and HIV) with high contamination:

0.5% 15 min.

Rotaviruses (EN 14476, clean conditions):

**0.5% 15 min.**

Noroviruses (EN 14476, high contamination):

1.5% 15 min.

1.0% 60 min.

Polyoma SV 40 viruses (EN 14476, dirty conditions):

1.0% 15 min.

#### **Extra details about the scope of effectiveness:**

- Bactericidal (incl. MRSA) as per EN 13727, EN 13697, EN 16615, VAH; levurocidal (*Candida albicans*) as per EN 13624, EN 13697, EN 16615, VAH; tuberculocidal as per EN 14348; limit virucidal effect (enveloped viruses such as HBV, HIV and HCV) as per DVV/RKI; effective against rotaviruses and noroviruses, as per EN 14476.

## **Safety guidelines**

H 302 Harmful if swallowed. H 314 Causes severe skin burns and eye damage. H 410 Very toxic to aquatic life with long lasting effects. Wear protective gloves/eye protection. Avoid release to the environment. IN CASE OF CONTACT WITH SKIN (or hair): Remove all contaminated clothing immediately. Wash/shower the skin with water. IN THE EVENT OF CONTACT WITH THE EYES: Rinse carefully for a long time with water. Remove any contact lenses if possible.

## **Available containers**

- Practical dispenser bottles
- Canisters

## **Contact**

PRISMAN Pharma International AG  
Am Stalden 16  
4622 Egerkingen  
Switzerland

Tel.: +41 (0)62 206 20 60

Fax: +41 (0)62 206 20 69

[www.innocid.com](http://www.innocid.com)

[info@prismanpharma.com](mailto:info@prismanpharma.com)

## **Manufacturer**

PRISMAN GmbH  
Otto-Hahn-Ring 6-18  
64653 Lorsch  
Germany

Continue to rinse. Call a POISON CENTER or a doctor/physician immediately. Dispose of the contents/container in accordance with local/regional/national/international regulations. Contains didecyldimethylammonium chloride, N-(3-aminopropyl)-N-dodecylpropane-1.3-diamin.



**Danger**