



# Product information

## Instrument Disinfection PREMIUM HIGH conc.

1000 ml	SUN404I
2500 ml	SUN404K
5000 ml	SUN404L



## Product type

For the disinfection of surgical instruments, as well as endoscopes of all kind.

## Application

The concentrate has to be diluted with water before use. According to the recommendation of RKI/BfArM, a thorough pre-cleaning or mechanical cleaning, as well as any necessary interim rinsing of the instruments before disinfection is to be ensured. Immerse the cleaned instruments into the (diluted) solution in an open condition. All surfaces and cavities to be disinfected have to be covered in full by the solution. Cover the baths used for the instrument disinfection. After the recommended application time, remove the instruments, thoroughly rinse with water and then dry, a final rinse with demineralised water is recommended. Disinfection solutions used are to be replaced at least every working day. When first used, or after a change of product, the instruments, baths and application aids are to be thoroughly cleaned and rinsed. Do not mix with other products (e.g. anionic surfactants, amines).

## Application times

Bacteria as per EN 13727, EN 14561 and VAH, with dirty conditions:

1.5% 15 min.

1.0% 60 min.

Candida albicans as per EN 13624, EN 14562, and VAH, with dirty conditions:

1.5% 15 min.

1.0% 60 min.

TbB (tuberculosis) as per EN 14348 and EN 14563 (dirty conditions):

3.0% 15 min.

1.5% 60 min.

Enveloped viruses (e.g. HBV, HCV and HIV) as per RKI/DVV (high contamination):

## Chemical-physical data

Density: 1.0 g/ml

pH-value: 4.5

## Available containers

- Practical dispenser bottles with dispenser container
- Canisters

0.5% 5 min.

Non-enveloped viruses (adenoviruses, polioviruses and noroviruses) as per EN 14476 (dirty conditions):

3.0% 15 min.

1.5% 60 min.

### Safety guidelines:

H 302 Harmful if swallowed. H 315 Causes skin irritation. H 318 Causes severe eye damage. H 334 May cause allergies, asthma symptoms or breathing difficulties if inhaled. H 317 May cause allergic skin reactions. H 335-H 336 May cause respiratory irritation. Can cause drowsiness or dizziness. H 400 Very toxic to aquatic life. Wear protective gloves/eye protection. IN THE EVENT OF CONTACT WITH THE EYES: Rinse carefully for a long time with water. Remove any contact lenses if possible. Continue to rinse. Call a POISON CENTER or a doctor/physician immediately. In case of symptoms affecting the breathing: Call a POISON CENTER or a doctor/physician. IF SWALLOWED: If you feel unwell contact the TOXIC INFORMATION CENTRE or a doctor/physician. Avoid release to the environment. Dispose of the contents/container in accordance with local/regional/national/international regulations. Contains glutaral (ethylenedioxy)dimethanol



**Danger**

### Certificates

Our products have been scientifically examined and tested in accordance with the latest methods. We will be happy to send you copies of our assessments on request.

### Contact

PRISMAN Pharma International AG

Am Stalden 16

4622 Egerkingen

Switzerland

Tel.: +41 (0)62 206 20 60

Fax: +41 (0)62 206 20 69

[www.innocid.com](http://www.innocid.com)

[info@prismanpharma.com](mailto:info@prismanpharma.com)

### Manufacturer

PRISMAN GmbH

Otto-Hahn-Ring 6-18

64653 Lorsch

Germany