



Product information

Surface Disinfection PREMIUM conc.

1000 ml	SUN310I
5000 ml	SUN310L
10000 ml	SUN310N



Characteristics

- wide microbiological spectrum of effectiveness
- no quaternary ammonium salts, particularly suitable for food purposes
- good material compatibility
- also tested for low temperatures (DVG)

Areas of application

Disinfectant concentrate for cleaning and disinfecting surfaces in the food processing industries as well as in the medical and institutional fields.

Instructions for use

Dilute the concentrate with tap water. Do not mix with cleaning agents or other disinfectants! Before using the product on large areas of sensitive materials like artificial leather, first test the material tolerance on a hidden area.

Wipe the treated areas at regular intervals with a damp disposable cloth (with water).

Application times

Bacteria and Candida Albicans (EN 16615, VAH, clean conditions):

1.0% 15 min.

Bacteria (dirty conditions):

1.0% 5 min. (EN 1276)

1.5% 15 min. (EN 13697)

1.0% 5 min. (EN 13727)

Bacteria as per DVG (8th list) at 20°C:

0.75% 5 min. (clean conditions)

0.25% 30 min. (clean conditions)

Composition

100g of the impregnating fluid contains:

6.75g N-(aminopropyl)-N-dodecylpropane-1.3-diamin.

Chemical-physical data

Density: 1.03g/ml

pH value: 10

1.0% 5 min (dirty conditions)

0.5% 30 min (dirty conditions)

Candida albicans (dirty conditions):

0.5% 15 min. (EN 1650)

1.0% 15 min. (EN 13697)

1.0% 5 min. (EN 13624)

Candida albicans as per DVG (8th list) at 20°C:

0.75% 30 min. (clean conditions)

1.0% 30 min. (dirty conditions)

Mycobacteria (*M. terrae* and *M. avium*) as per EN 14348:

1.5% 60 min. (clean conditions)

2.0% 60 min. (dirty conditions)

Enveloped viruses (e.g. HBV, HCV and HIV) as per EN 14476, high contamination:

2.0% 15 min.

1.5% 30 min.

Safety instructions

H290 May be corrosive to metals. H 315 Causes skin irritation, H 318 Causes serious eye damage. H412

Very toxic to aquatic life with long-lasting effects.

Wear protective gloves/eye protection/face

protection. Avoid release to the environment. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor/... Absorb spillage to prevent material damage. Dispose of contents/container into recognised waste disposal system.

Contains N-(3-aminopropyl)-N-dodecylpropane-1.3-diamin.

Only put the packaging into the recycling system if it is empty of residue. Dispose of residual product as per EAK 07 06 01.

Delivery options

- Practical dispenser bottles with dispenser container
- Canisters

Certificates

Our products have been scientifically examined and tested in accordance with the latest methods. We will be happy to send you copies of our assessments upon request.

Contact

PRISMAN Pharma International AG
Am Stalden 16
4622 Egerkingen
Switzerland

Tel.: +41 (0)62 206 20 60

Fax: +41 (0)62 206 20 69

www.sunsept.de

info@prismanpharma.com

Manufacturer

PRISMAN GmbH
Otto-Hahn-Ring 6-18
64653 Lorsch
Germany



Use biocidal products cautiously. Always read label and product information before use.