



Product information

Surface Disinfection PLUS Konz.

1000 ml	SUN370I
2500 ml	SUN370K
5000 ml	SUN370L
10000 ml	SUN370N



- Excellent cleaning performance (no additional cleaning additives necessary)
- Does not leave sticky residues
- Extremely broad spectrum of efficacy with low dosage
- Also suitable for dosing systems
- Aldehyde- and phenol-free

Product-type

Water-based disinfectant concentrate for cleaning and disinfecting surfaces. A comprehensive spectrum of activity against bacteria, yeasts and viruses is covered at low application concentrations.

Areas of application

As a biocidal product:

Aldehyde-free surface disinfection concentrate for cleaning and disinfection of surfaces.

As a medical device:

Surface disinfection for the prophylactic surface disinfection of medical devices. Not suitable for the final disinfection of invasive medical devices.

Instructions for use

Before use, dilute the disinfectant concentrate with tap water. Do not mix with cleaning agents or other disinfectants! Before large-scale applications on sensitive materials such as artificial leather, test for material compatibility on a small hidden area. To avoid streaks, wipe the treated surfaces at regular intervals with a damp (water) disposable wipe. Use only as wipe disinfection! For professional users only. Fresh ready-to use solution should be made up daily.

Composition

100 g contain:

10.8 g Dodecyl dimethyl ammonium chloride

6.0 g N-(3-aminopropyl)-N-dodecylpropan-1,3-diamine

Chemical-physical data

Density: 1.01 g/m³

pH value (at 20° C) >11

Application time

Bacteria and Candida albicans (EN 13727, EN 13697, EN 13624, VAH, EN 16615 in dirty conditions):

0.5 % 15min.

TB (*M. terrae*) and EN 14348:

2.0 % 15min. (clean conditions)

0.5 % 60min. (clean conditions)

1.0 % 60min. (dirty conditions)

„Limited virucidal“ (enveloped viruses as per RKI/DVV such as HBV, HCV, HIV) dirty conditions:

0.5 % 15min.

Rotaviruses (EN 14476, clean conditions):

0.5 % 15min.

Noroviruses (EN 14476, dirty conditions):

1.5 % 15min.

1.0% 60min.

Adenoviruses (ENB 14476, clean conditions):

2.0 % 30 min.

“Limited virucidal PLUS” (EN 14476, clean condition):

2.0 % 30 min.

Polyomavirus (SV40) (EN 14476, dirty conditions):

1.0 % 15min.

Additional details on the efficacy spectrum

Bactericidal (incl. MRSA) as per EN 13727, EN 13697, EN 16615, VAH; yeasticidal (*Candida Albicans*) as per EN 13624, EN 13697, EN 16615, VAH; tuberculocidal as per EN 14348; “limited virucidal Plus” (enveloped viruses such as HBV, HIV, HCV) as per DVV/RKI and noro-, adeno- and rotaviruses as per EN 14476).

Available containers

- Practical dispenser bottles
- Canisters

Contact

PRISMAN Pharma International AG
Am Stalden 16
4622 Egerkingen
Switzerland

Tel.: +41 (0)62 206 20 60

Fax: +41 (0)62 206 20 69

www.innocid.com

info@prismanpharma.com

Manufacturer

PRISMAN GmbH
Otto-Hahn-Ring 6-18
64653 Lorsch
Germany

Safety instructions

H 302 Harmful if swallowed. H 314 causes severe skin burns and eye damage. H 410 Very toxic to aquatic life with long lasting effects. Wear protective gloves / eye protection. Avoid release to the environment. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor. Dispose of contents/container in accordance with local/regional/national/international regulations. Contains didecyl dimethyl ammonium chloride, N-(3-Aminopropyl)-N-dodecyl-propan-1,3-diamine

